

PHOSPHORUS (BOUND)

PRINCIPLE

The sample is washed thoroughly with an aqueous solvent to remove soluble phosphorus-containing materials, and then ignited in the presence of a fixative to destroy organic matter and convert phosphorus to inorganic phosphates which are not volatilized during ignition. Residual phosphates are taken up in acid, hydrolyzed to orthophosphate and determined spectrophotometrically as the molybdovanadophosphoric acid complex (Note 1).

SCOPE

The procedure is applicable to the determination of bound phosphorus in unmodified, modified and derivatized granular starches. With only minor modification, the procedure can be applied to other cold water insoluble products obtained from the wet milling of corn and grain sorghum (Note 2).

SPECIAL APPARATUS

1. VYCOR Dishes: 100 mL capacity - Do not use platinum (Note 3).
2. Muffle Furnace: Equipped with a pyrometer and capable of operating at controlled temperatures up to 650 °C
3. Double Beam Spectrophotometer: An instrument capable of accurate absorbance measurements at 460 nm and equipped with a matched pair of 1 cm cuvettes.

REAGENTS

1. Methanol-Water, 7:3 (V/V): Add 700 mL of reagent grade methanol to 300 mL of purified water and mix.
2. Zinc Acetate Solution, 10%: Dissolve 120.0 g of reagent grade zinc acetate dihydrate ($\text{Zn}(\text{C}_2\text{H}_3\text{O}_2)_2 \cdot 2\text{H}_2\text{O}$) in 880 mL of purified water (Note 4). If the solution is hazy, filter through Whatman No. 2V filter paper.

**Analytical Methods of the Member Companies of the
Corn Refiners Association, Inc.**

Accepted 10-16-69
Revised 10-8-96

PHOSPHORUS (BOUND) — continued

3. Standard Phosphorus Solution, 0.1 mg/mL as Phosphorus: Dissolve exactly 0.4395 g (Note 5) of reagent grade potassium dihydrogen phosphate (KH_2PO_4) in water and dilute to 1 L.
4. Ammonium Vanadate Solution, 0.25%: Dissolve 2.5 g of ammonium metavanadate ($(\text{NH}_4)\text{VO}_3$) in 600 mL of boiling purified water. Cool to 60-70 °C and add 20 mL of concentrated nitric acid. Cool to room temperature, transfer quantitatively to a 1 L volumetric flask and dilute to mark with purified water.
5. Ammonium Molybdate Solution, 5%: Dissolve 50 g of ammonium molybdate tetrahydrate ($(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot\text{H}_2\text{O}$) in 900 mL of warm purified water. Cool to room temperature and dilute to 1 L with purified water.
6. Nitric Acid Solution, 29%: Add 300 mL of concentrated nitric acid (70% HNO_3) to 600 mL of purified water and mix.

PROCEDURE

Standardization: Pipet 5.0, 10.0, 15.0, 20.0 and 25.0 mL of standard phosphorus solution (respectively, 0.5, 1.0, 1.5, 2.0 and 2.5 mg of phosphorus) into respective 100 mL volumetric flasks and use another flask for a blank. To each flask add, *in order*, 10 mL of 29% nitric acid, 10.0 mL of 0.25% ammonium vanadate, and 10 mL of 5% ammonium molybdate, mixing thoroughly after addition of each reagent (Note 6). Dilute to volume with water, mix thoroughly and allow to stand for 10 mins. (Note 7). Using the blank as a reference solution at zero absorbance, determine absorbance of each standard at 460 nm. Plot absorbance vs. mg of phosphorus per 100 mL (Note 8).

Sample Preparation: Place 20 to 25 g of granular starch product (Note 9) in a 250 mL beaker and add 200 mL of 7:3 methanol-water at room temperature. Disperse starch and agitate mechanically for 15 mins. Recover starch by vacuum filtration in a 150 mL medium porosity fritted glass or Büchner funnel and wash the wet cake with 200 mL of 7:3 methanol-water. Reslurry and rewash the wet cake with 7:3 methanol-water.

PHOSPHORUS (BOUND) — continued

Dry the washed starch filter cake in an air oven at a temperature below 50 °C. Grind the sample completely through a laboratory cutting mill or with mortar and pestle to 20 mesh or finer and blend thoroughly. Determine dry substance by a standard procedure.

Low-Phosphorus Samples (0.005% to 0.05% P): Weigh accurately in a VYCOR dish an amount of washed and dried starch sample containing not more than 2.5 mg of phosphorus (Note 3). Recommended sample sizes for low-phosphorus starches and other sample types are approximately as follows:

<u>Phosphorus, % as is</u>	<u>Sample Wt. (g)</u>
0.01 and lower	20
0.01 to 0.025	10
0.025 to 0.050	5

High-Phosphorus Samples (0.05% to 2.0% P): Weigh accurately an amount of washed and dried starch sample containing not more than 40 mg of phosphorus (Note 9). Recommended sample sizes for high phosphorus starches and other sample types are listed below, together with recommended residue dilution volumes and aliquot sizes for spectrophotometric analysis.

<u>Phosphorus, % as is</u>	<u>Sample Wt. (g)</u>	<u>Dilution Volume (mL)</u>	<u>Aliquot Volume (mL)</u>
0.05 to 0.10	10	100	25
0.10 to 0.25	10	200	20
0.25 to 0.75	5	250	15
0.75 to 2.00	2	250	15

In all cases, add 10 mL of 10% zinc acetate solution to the VYCOR and distribute the solution uniformly through the sample adding water if necessary. Evaporate to dryness on a steam bath and char the sample on a hot plate. Ignite in a muffle furnace at 550 °C for 2 hrs.

PHOSPHORUS (BOUND) — continued

Cool to room temperature and wet all of the residue by cautious addition of 3 mL of 29% nitric acid (Note 10). Evaporate to dryness on a steam bath and dehydrate by brief heating on a hot plate. Return the dish to the muffle furnace for 20 mins.

Cool to room temperature, wash down the sides of the dish with 10 mL of 29% nitric acid and add 15 mL of purified water. Cover with a watch glass, heat to incipient boiling and hold at that temperature for 10 mins. (Note 11).

For low phosphorus samples: Quantitatively filter through Whatman No. 1 filter paper into a 100 mL volumetric flask. Complete the transfer with the aid of four 10 mL portions of purified water. To each flask, add *in order*, 10 mL of 29% nitric acid, 10.0 mL of 0.25% ammonium vanadate, and 10 mL of 5% ammonium molybdate, mixing thoroughly after addition of each reagent (Note 6).

For high phosphorus samples: Filter into a volumetric flask of a size selected to yield a final phosphorus concentration of 1 to 4 mg/mL P. Complete the transfer with the aid of four 10 mL portions of purified water. Dilute to volume with purified water and mix thoroughly. Transfer an aliquot selected to contain not more than 2.5 mg of phosphorus to a 100 mL volumetric flask and add 50 mL of purified water to another flask to serve as a blank. Add, *in order*, 10 mL of 29% nitric acid, 10 mL of 0.25% ammonium vanadate, and 10 mL of 5% ammonium molybdate, mixing thoroughly after addition of each reagent (Note 6).

In all cases, dilute to volume with purified water, mix thoroughly and allow to stand for 10 mins. (Note 7). Determine absorbance of the sample at 460 nm using the blank as a reference solution at zero absorbance. Read mg of phosphorus in the aliquot from the standardization curve.

CALCULATIONS

Low-Phosphorus Samples:

$$\% \text{ (dry basis) Phosphorus} = \frac{(\text{Phosphorus, mg from Graph})(100)}{(\text{Dry Sample Wt., g})(1000)}$$

PHOSPHORUS (BOUND) — continued

High-Phosphorus Samples:

$$\% \text{ (dry basis)} = \frac{(\text{Phosphorus, mg from Graph})(\text{Dilution Volume})(100)}{(\text{Aliquot Volume})(\text{Dry Sample Wt., g})(100)}$$

NOTES AND PRECAUTIONS

1. Chemistry of the chromophore is in dispute; for a discussion of some of the possibilities, see C.Y. Shen and D.R. Dyroff, Determination of Phosphate in Presence of Silicates by Molybdenum Blue Method; *Anal. Chem.*, **1962**, Vol. 34, no. 11, pp. 1367-1370.
2. The procedure as written is not recommended for starches containing less than 40 ppm phosphorus. In this low range, adequate precision can be attained by ten-fold dilution of the phosphorus standard, and use of 10 cm cuvettes for analysis of standards and samples.
3. Platinum *must not* be used, because it may be damaged by the fixative or by large phosphate residues.
4. Zinc acetate dihydrate has not been observed to contribute apparent phosphorus to the system; however, it is advisable to check an ignited zinc acetate blank against a reagent blank whenever a new lot of zinc acetate dihydrate is opened.
5. Most samples of potassium dihydrogen phosphate declared assay slightly over 100% KH_2PO_4 . Obtain the exact weight of standard to be taken by dividing 0.4395 by the decimal purity factor.

(Example: declared assay = 100.4% KH_2PO_4 ; decimal purity factor = 1.004; weight of standard = 0.4377 g)
6. Reagents *must* be added *in order* as stated in order to avoid interferences from precipitation and side reactions.
7. The color increases very slowly for the first four hours, resulting in a relative error of about 2%. With double beam spectrophotometers, error is

PHOSPHORUS (BOUND) — continued

negligible if readings are made within 2 hrs. of solution preparation or if samples and standards are read at comparable intervals after preparation.

8. The standardization curve is reproducible and needs be checked only when fresh reagents are prepared.
9. If the sample is soluble in cold water, prepare a 1 to 2% aqueous solution (paste); place in a cellophane tube and dialyze against running purified water for 30-40 hrs. Precipitate the starch by pouring the solution into 4 volumes of acetone (per volume of paste) while stirring. Recover precipitated product by vacuum filtration on a fritted glass or Büchner funnel. Wash filter cake with absolute ethanol. Dry in an air or vacuum oven, grind, blend and determine volatiles before analysis. A dry product can also be recovered from the dialyzed paste by freeze-drying.
10. The residue will seldom be carbon-free at this point; treatment with nitric acid and re-ignition assures complete combustion in a reasonable time. Residue from the evaporation contains hydrates which must be decomposed, by heating as directed, to avoid loss by spattering in the furnace.
11. Part or all of the phosphorus in the residue is present as pyrophosphate, which will not form the molybdovanadophosphoric acid complex. The specified heating period is necessary to hydrolyze pyrophosphate to the reactive orthophosphate. Since higher phosphates, such as trimetaphosphate or tripolyphosphate, may be present (as in unreacted mixtures of starch with polymeric phosphates), it is advisable to extend the heating time to 30 mins.